

Integrating Device Engineering and User-Centered Design for High-Performance Dry Powder Inhalers Using a Particle-Engineered Formulation

Ângelo Araújo¹, Sofia Silva¹, Filipe Vultos¹, Mirjam Kobler², Ameet Sule², Sunita Sule², David Valencia Pellisa², Daniel Garcia Sabido²
¹R&D, Hovione FarmaCiência SA, Lumiar, 1649-038, Portugal. ²Presspart Manufacturing Limited, Blackburn, BB15RF.

BACKGROUND

The Sunriser[®] DPI was developed using an integrated approach and outperformed a leading marketed DPI in both technical and ergonomic metrics




Figure 1. Sunriser[®] DPI

- DPIs are popular for respiratory drug delivery due to **convenience and stability** but **achieving consistent aerosol performance** while keeping devices **user-friendly** is challenging.
- Carrier-based blends** work at **low doses** but fail at high loads due to weak interactions and cohesion; **carrier-free spray-dried powders** suit **high doses** but often agglomerate and aerosolize poorly, especially for biologics and vaccines and emergency therapies.
- A stepwise, feedback-driven development strategy—using a spray-dried for **device engineering**, with rapid prototyping (e.g., 3D printing)—enables iterative testing and refinement, with real-time usability and lab feedback **improving both technical performance and user experience**.

The objective of this work is to develop a high-performance, user-friendly DPI by integrating formulation science, device engineering, and user feedback to overcome current limitations and enhance lung deposition and patient compliance.


METHODS



RESULTS AND DISCUSSION

USER-CENTERED DESIGN

- 30-participant survey assessed ergonomics (grip, cap, button placement).
- Feedback informed refinements in Prototype III:



PROTOTYPE EVOLUTION





Figure 2. FPD and ED of the marketed device compared with the three-stage prototypes. Asterisks indicate significant pairwise differences (*p < 0.05; **p < 0.01; ***p < 0.001).

Table 1. Analytical performance of three prototype stages using an engineered formulation, compared with marketed device (LC = label claim, N=3)

	Marketed device	Prototype I Initial 3D printed	Prototype II Optimized 3D printed	Prototype III Representative
ED (% LC ± SD)	80.4 ± 0.45	72.9 ± 2.12	78.6 ± 1.31	80.5 ± 4.32
FPD (mg ± SD)	10.5 ± 0.88	12.5 ± 0.48	13.7 ± 0.14	13.6 ± 0.50

CONCLUSIONS



PROTOTYPE EVOLUTION

The Sunriser[®] DPI demonstrates how a structured, multidisciplinary development strategy can overcome limitations of conventional DPIs.



FORMULATION SCIENCE, DEVICE ENGINEERING, AND USER FEEDBACK

Integration of formulation science, device engineering, and user feedback led to a production-ready inhaler with superior performance and usability.



FUTURE DPI INNOVATION

This approach sets a benchmark for future DPI innovation, especially for high-dose and biologic therapies.

References

- Groß R, et al., J Pharmaceutics 2022; 14 pp1185–1199.
- Omidian H, et al., Pharmaceuticals 2025; 18(2); pp 175.
- Ye, Y, et al., Int J Pharmaceutics 2022; 25:614:121457
- Saldanha S, et al., DDL 2023.
- Abrego Morales L, et al., Pharmaceutics 2025; 17(3); pp 359.
- Kobler M, et al., ISAM 2025.
- Churro R, et al., RDD 2025; pp 75–79.
- Kobler M, et al., Respiratory Drug Delivery 2024; 1: 231-234.